


IFC PAL Fraction Collection at Analytical and Prep Scale




Collect Sample


Empty Syringe

Divert to Waste

Sideport Syringe Function



IFC PAL equipped for analytical applications



IFC PAL equipped for preparative applications

- Sample injection and collection on one single instrument
- Collects in up to 1000 test tubes or 24 microplates
- Dual arm concept prevents cross contamination risk
- Unique side hole syringe acts as diverter valve
- Variable rack layout allows quick process adaption
- Optimized time management, next injection prepared during collection
- One software for injection, collection and LC pumps